

# CLINICAL OUTCOME DATA FOLLOWING ROTATOR CUFF REPAIR USING SINEFIX™

## A RETROSPECTIVE ANALYSIS OF FUNCTIONAL AND PAIN-RELATED PARAMETERS

Stefan Welte, MD, Tobias Goetz, Keith Hall, MD

### EXECUTIVE SUMMARY

This paper summarizes retrospective clinical outcome data from patients undergoing rotator cuff repair using SINEFIX, with a focus on postoperative pain development during early follow-up.

### METHODS

A retrospective review was conducted in 20 consecutive patients who underwent rotator cuff repair using SINEFIX. Pain intensity was assessed using a visual analogue scale (VAS) preoperatively and during routine postoperative follow-up visits up to approximately 90 days. n = 19 patients had complete pre- and postoperative VAS data available for analysis.

### RESULTS


In the analyzed cohort (n = 19), mean VAS pain scores decreased significantly from 5.89 (SD 2.49) preoperatively to 2.53 (SD 2.65) at the last postoperative follow-up (Wilcoxon signed-rank test,  $p < 0.01$ ), indicating a clinically meaningful reduction in patient-reported pain during the early postoperative period.

No intraoperative complications were reported. Furthermore, no patient expressed dissatisfaction with their outcome during the available follow-up period. The operating surgeon also reported that the new technique was straightforward to perform and easy to integrate into the surgical workflow.

### CONCLUSION

In this retrospective dataset, rotator cuff repair using SINEFIX was associated with reductions in pain within the studied patient cohort.

No intraoperative complications were observed, and patient satisfaction was consistently positive. Further prospective studies including functional and structural outcome measures are warranted to confirm these findings.



**FIGURE 1**

Mean preoperative versus last follow-up VAS pain scores following rotator cuff repair using SINEFIX. Retrospective analysis, n = 19.

---

**Inovedis Inc.**

7625 Golden Triangle Dr., Suite G  
Eden Prairie, MN 55344

Phone: 612-445-8362

Mail: [cs@inovedis.com](mailto:cs@inovedis.com)

Legal Manufacturer: BAAT Medical Products B.V., F. Hazemeijerstraat 800 - A04 | 7555 RJ Hengelo | The Netherlands

**DISCLAIMER**

This white paper is for informational purposes only and does not constitute medical, legal, or professional advice. The information contained herein is based on the knowledge and research available at the time of writing and is subject to change without notice. The authors and publishers of this white paper make no representations or warranties regarding the completeness, accuracy, reliability, or suitability of the information contained in this document. Readers are advised to consult with qualified medical professionals for advice, diagnosis, or treatment related to implants or any other medical conditions. The content of this white paper should not be used as a substitute for professional medical consultation or treatment. Any reliance on the information provided in this document is solely at the reader's own risk. The authors, publishers, and affiliated entities disclaim any liability for any direct, indirect, incidental, or consequential loss or damage incurred by individuals or entities relying on the information presented in this white paper. By using this white paper, the reader acknowledges and agrees to the terms of this disclaimer.

SINEFIX is FDA cleared and approved for use in the U.S.

MKT-WP-20260127 Rev B